

50A7 Epi Sensa

DE	Gebrauchsanweisung	3
EN	Instructions for use	7
FR	Instructions d'utilisation	11
IT	Istruzioni per l'uso	15
ES	Instrucciones de uso	19
PT	Manual de utilização	23
NL	Gebruiksaanwijzing	28
SV	Bruksanvisning	32
DA	Brugsanvisning	35
NO	Bruksanvisning	39
FI	Käyttöohje	43
PL	Instrukcja użytkowania	47
HU	Használati utasítás	51
CS	Návod k použití	55
HR	Upute za uporabu	59
SK	Návod na používanie	63
BG	Инструкция за употреба	67
TR	Kullanma talimatı	71
EL	Οδηγίες χρήσης	75
RU	Руководство по применению	79
ZH	使用说明书	84

	Größe Size	Unterarmumfang Forearm circumference	
		cm	inch
	XS	17 - 19	6.7 - 7.5
	S	19 - 21	7.5 - 8.3
	M	21 - 23	8.3 - 9.0
	L	23 - 25,5	9.0 - 10.0
	XL	25,5 - 28	10.0 - 11.0
	XXL	28 - 31	11.0 - 12.2

Material	Gestrick / Knit Fabric: PA, EL
	Frotteebezüge / Terry Cloth Covers: Baumwolle / Cotton, PES, EL
	Pelotte / Pad: viskoelastisches Technogel® / viscoelastic Technogel®

INFORMATION

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INFORMATION

Date of last update: 2022-12-20

- ▶ Please read this document carefully before using the product and observe the safety notices.
- ▶ Instruct the user in the safe use of the product.
- ▶ Please contact the manufacturer if you have questions about the product or in case of problems.
- ▶ Report each serious incident related to the product to the manufacturer and to the relevant authority in your country. This is particularly important when there is a decline in the health state.
- ▶ Please keep this document for your records.

These instructions for use provide important information on the fitting and application of the 50A7 Epi Sensa elbow support.

2 Intended use

2.1 Indications for use

The support is intended **exclusively** for treating the upper limbs and **exclusively** for contact with intact skin.

The support must be used according to the indication.

2.2 Indications

- Epicondylitis humeri radialis (tennis elbow)
- Epicondylitis humeri ulnaris (golfer elbow)
- Tendomyopathies and other soft tissue disorders (e.g. chronic bursitis) and injuries
- Joint effusion and periarticular swelling (e.g. in case of osteoarthritis and arthritis)
- Postoperative and posttraumatic irritation

The indication must be determined by the physician.

2.3 Contraindications

2.3.1 Absolute Contraindications

None known.

2.3.2 Relative Contraindications

The following indications require consultation with a physician: skin diseases/injuries, inflammation, prominent, swollen scars, reddening and hyperthermia in the relevant area; pronounced varicose veins, especially with impaired blood return, lymph drainage disorders and unclear soft tissue swelling distal to the device; sensory and circulatory disorders in the area of the forearm and hand, e.g. associated with diabetic neuropathy.

2.4 Lifetime

The product is designed for a maximum lifetime of **6 months**.

2.5 Mechanism of Action

The form-knitted support with integrated shape-retaining pad provides elbow compression. The support promotes sensorimotor functions and thereby activates muscular stabilization.

3 Safety

3.1 Explanation of warning symbols

CAUTION Warning regarding possible risks of accident or injury.

NOTICE Warning regarding possible technical damage.

3.2 General safety instructions

Contact with heat, embers or fire

Risk of injury (such as burns) and risk of product damage

- Keep the product away from open flames, embers and other sources of heat.

Reuse on other persons and improper cleaning

Skin irritation, formation of eczema or infections due to contamination with germs

- The product may be used by one person only.
- Clean the product regularly.

NOTICE

Contact with oils, salves, lotions or other products that contain oils or acids

Insufficient stabilization due to loss of material functionality

- ▶ Do not expose the product to oils, salves, lotions or other products that contain oils or acids.

4 Handling

INFORMATION

- ▶ The daily duration of use and period of application are generally determined by the physician.
- ▶ The initial fitting and application of the product must be carried out by qualified personnel.
- ▶ Instruct the patient in the handling and care of the product.
- ▶ Instruct the patient to see a physician immediately if any exceptional changes are noted (e.g. worsening of the complaint).

4.1 Selecting the Size

- 1) Measure the circumference of the forearm (approximately 12 cm below the elbow).
- 2) Determine the size of the support (see size chart).

4.2 Application

⚠ CAUTION

Incorrect or excessively tight application

Skin irritation due to overheating, local pressure points due to tight application to the body

- ▶ Ensure that the product is applied properly and fits correctly.
- ▶ When in doubt, do not continue to use the product if skin irritation occurs.

NOTICE

Use of a worn or damaged product

Limited effectiveness

- ▶ Before each use, check the product for functional reliability and for possible wear or damage.

- Do not continue using a product that is no longer functional, or that is worn or damaged.

- 1) Fold over the top section of the support so that the pad pockets are exposed (see fig. 1).
- 2) Pull the support over the forearm up to the elbow.
- 3) Fold back the section of the support that was folded over and pull the support over the elbow (see fig. 2).
- 4) Position the support so that the epicondylus lateralis and epicondylus medialis are positioned in the pad openings (see fig. 3).
- 5) When the patient is given the support, it must be checked for proper fit and positioning.

5 Cleaning

NOTICE

Use of improper cleaning agents

Damage to the product due to use of improper cleaning agents

- Only clean the product with the approved cleaning agents.

Clean the support regularly:

- 1) **Recommendation:** Use a laundry bag or net.
- 2) Wash the support in warm water at 30 °C (86°F) with standard mild detergent. Do not use fabric softener. Rinse thoroughly.
- 3) Allow to air dry. Do not expose to direct heat (e.g. sunshine, stove or radiator).

6 Disposal

Dispose of the product in accordance with national regulations.

7 Legal Information

All legal conditions are subject to the respective national law of the country of use and may vary accordingly.

7.1 Liability

The manufacturer shall be liable in the event that the product is used in accordance with the descriptions and instructions in this document. The manufacturer will not assume liability for damage caused by disregarding the information in this document, particularly due to improper use or unauthorized modification of the product.

7.2 CE Conformity

The product meets the requirements of Regulation (EU) 2017/745 on medical devices. The CE declaration of conformity can be downloaded from the manufacturer's website.

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